




DM1&DM4 Recombinant Monoclonal Antibody

Product Code	CSB-RA989811I2m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	DM4-BSA
Species Reactivity	Chemical
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG2a
Clonality	Monoclonal
Product Type	Recombinant Antibody
Research Area	Others
Target Names	DM1&DM4
Clone No.	12H10E10

Image



The Binding Activity of Trastuzumab emtansine (Ado-Trastuzumab emtansine) and Cantuzumab ravtansine (Synonyms:IMGN242) with Anti-DM1&DM4 antibody

Activity: Measured by its binding ability in a functional ELISA. Immobilized Trastuzumab emtansine (Ado-Trastuzumab emtansine) and Cantuzumab ravtansine (Synonyms:IMGN242) at 2 µg/mL can bind Anti-DM1&DM4 antibody, the EC₅₀ is 6.330 to 8.807 ng/mL and 10.49 to 15.04 ng/mL, respectively.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.