




Recombinant Mouse Cytochrome P450 2C55 (Cyp2c55), partial

Product Code	CSB-EP875232MO
Abbreviation	Recombinant Mouse Cyp2c55 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9D816
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Mouse Cytochrome P450 2C55 (Cyp2c55), partial
Immunogen Species	Mus musculus (Mouse)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MDPVLVLVLTLSCLLLLSLWRQNSGRGKLPPGPTPFPIIGNILQIDIKNISKSFNY FSKVYGPVFTLYFGSKPTVVVHGYEAVKEALDDLGEETFSGRGSFQIFERINNDL GVIFSNGTKWKELRRFSIMTLRSFGMGKRSIEDRIQEEASCLVEELRKANGSLC DPTFILSCAPSNVICSVIFHNRFDYKDEKFLNLMERLNENFKILNSPWMQVYNA LPTLINYLPGSHNKVIKNFTEIKSYILGRVKEHQETLMDNPRDFIDCFLIKMEQE KHNPHEFTIESLMATVTDIFVAGTETTITLRYGLLLLLKHTEVTAKVQAEIDHV IGRHRSPCMQDRTRMPYTDAMVHEIQRYIDLIPNNVPHAATCNVRFRSYFIPK GTELVTSLTSVLHDDKEFPNPEVFDPGHFLDENGNFKKSDYFMPFSIGKRMCV GEALARTEFLILTTILQNFNLKSLVDTKDIDTTPVANTFGRVPPSYQLYFIPR
Research Area	Cardiovascular
Source	E.coli
Target Names	Cyp2c55
Expression Region	1-490aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	63.0 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.