




# Recombinant Epstein-Barr virus Latent membrane protein 1 (LMP1)-VLPs

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-MP355987EFA                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Abbreviation</b>      | Recombinant Epstein-Barr virus LMP1 protein-VLPs                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                    |
| <b>Uniprot No.</b>       | P03230                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Form</b>              | Lyophilized powder                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Storage Buffer</b>    | Lyophilized from a 0.2 μm filtered PBS, 6% Trehalose, pH 7.4.                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Epstein-Barr virus Latent membrane protein 1 (LMP1)-VLPs                                                                                                                                                                                                                                                                                                                                                               |
| <b>Sensitivity</b>       | Not Test                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Purity</b>            | ≥ 90% as determined by SEC-HPLC.                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Sequence</b>          | MEHDLERGPPGPRRPPRGPPPLSSSLGLALLLLLLLALLFWLYIVMSDWTGGALL<br>VLYSFALMLIIIIILIFIFRRDLLCPLGALCILLMITLLIIALWNLHGQALFLGIVLFIF<br>GCLLVLGIWIYLLLEMLWRLGATIWQLLAFFLAFFDLILLIIALYLQQNWWTLLVD<br>LLWLLLFLAILIWMYYHGQRHSDEHHHDDSLPHPQQATDDSGHESDSNSNEG<br>RHHLLVSGAGDGPPLCSQNLGAPGGGPDNGPQDPDNTDDNGPQDPDNTDD<br>NGPHDPLPQDPDNTDDNGPQDPDNTDDNGPHDPLPHSPSDSAGNDGGPPQ<br>LTEEVENKGGDQGPPLMTDGGGGHSHDSGHGGGDPHLPPTLLLGS SSGSGGD<br>DDDPHGPVQLSYYD |
| <b>Research Area</b>     | Others                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Source</b>            | Mammalian cell                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Target Names</b>      | LMP1                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Expression Region</b> | 1-386aa                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Notes</b>             | The VLPs are expressed from human 293 cells (HEK293). Mix the sample gently by repeatedly pipetting it up and down. Do not vortex. Repeated freezing and thawing is not recommended. Store the protein at -20°C/-80°C upon receiving it, and ensure to avoid repeated freezing and thawing, otherwise, it will affect the protein activity. The immunization strategy should be optimized (antigen dose, regimen and adjuvant).    |
| <b>Tag Info</b>          | C-terminal 10xHis-tagged                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Mol. Weight</b>       | 43.4 kDa                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Protein Length</b>    | Full Length                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Image</b>             | <div>  </div> <div>CSB-MP355987EFA is detected by Mouse anti-6*His monoclonal antibody.</div>                                                                                                                                                                                                                                                   |



The purity of VLPs was greater than 90% as determined by SEC-HPLC

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. Aliquot for long-term storage at -80°C. Solubilize for 60 minutes at room temperature with occasional gentle mixing. Avoid vigorous shaking or vortexing.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.