




CD19 Recombinant Monoclonal Antibody

Product Code	CSB-RA004888MA4HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15391
Immunogen	Recombinant Human CD19 protein
Species Reactivity	Human
Tested Applications	FC; Recommended dilution: FC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology;Stem cells
Target Names	CD19
Clone No.	12G5

Image



Raji cells were stained with Human IgG1 Isotype Control (CSB-RA011156MA1HU)(green line) and anti-CD19 antibody (CSB-RA004888MA4HU) (2 μ g/1*10⁶cells) (red line), washed and then followed by FITC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.