





TSPAN8 Recombinant Monoclonal Antibody

Product Code	CSB-RA025166MA3HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19075
Immunogen	Recombinant Human TSPAN8 protein
Species Reactivity	Human
Tested Applications	FC; Recommended dilution: FC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	mIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 95% as determined by SEC-HPLC.
Research Area	Cancer;Tags & Cell Markers;Signal transduction
Target Names	TSPAN8
Clone No.	7D10

Image



Untransfected MC38 cells (green line) and transfected Human TSPAN8 MC38 stable cells (red line) were stained with anti-TSPAN8 antibody (CSB-RA025166MA3HU) ($2\mu\text{g}/1 \times 10^6$ cells), washed and then followed by FITC-conjugated anti-Mouse IgG Fc antibody and analyzed with flow cytometry.



The purity of TSPAN8 was greater than 95% as determined by SEC-HPLC



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.