





Recombinant Human Thrombopoietin (THPO) (Active)

Product Code	CSB-MP023509HU(A4)
Abbreviation	Recombinant Human THPO protein (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40225
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, 6% Trehalose, pH 7.4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Biological Activity	The ED50 as determined by the dose-dependent stimulation of the proliferation of MO7e cells is 0.4353 - 1.089 ng/mL.
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	SPAPPACDLRVLSKLLRDSHVLHSRLSQCPEVHPLPTPVLLPAVDFSLGEWKT QMEETKAQDILGAVTLLLEGVMAARGQLGPTCLSSLLGQLSGQVRRLLLGALQS LLGTQLPPQGRTTAHKDPNAIFLSFQHLLRGKVRFLMLVGGSTLCVRRAPPTT AVPSRTSLVLTNLNPNRTSGLLETNFTASARTTGSGLLKWQQGFRAKIPGLLN QTSRSLDQIPGYLNRIHELLNGTRGLFPGPSRRTLGAPEISSGTSDTGSLPPNL QPGYSPSPTHPPTGQYTLFPLPPTLPTPVVQLHPLLDPSPAPTPTPTSPLLNTS YTHSQNLSQEG
Source	Mammalian cell
Target Names	THPO
Expression Region	22-353aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	36.8 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Activity

The ED₅₀ as determined by the dose-dependent stimulation of the proliferation of MO7e cells is 0.4353 - 1.089 ng/mL.

Endotoxin

Less than 1.0 EU/ug as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.