





Recombinant Human B-cell antigen receptor complex-associated protein alpha chain (CD79A), partial (Active)

Product Code	CSB-MP004957HU1p6
Abbreviation	Recombinant Human CD79A protein, partial (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P11912
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, 6% Trehalose, pH 7.4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Biological Activity	Measured by its binding ability in a functional ELISA. Immobilized Human CD79A at 2 µg/mL can bind Anti-CD79A recombinant antibody(CSB-RA004957MA1HU). The EC50 is 14.44-16.80 ng/mL.
Purity	≥ 95% as determined by SDS-PAGE.
Sequence	LWMHKVPASLMVSLGEDAHFQCPSHNSNNANVTWWRVLHGNYTWPPEFLG PGEDPNGTLIIQNVNKS HGGIYVCRVQEGNESYQQSCGTYLRVRQPPRPFL DMGEGTKNR
Source	Mammalian cell
Target Names	CD79A
Expression Region	33-143aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal mFc-Flag-tagged
Mol. Weight	41.3 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Activity

Measured by its binding ability in a functional ELISA. Immobilized Human CD79A at 2 µg/ml can bind Anti-CD79A recombinant antibody (CSB-RA004957MA1HU). The EC₅₀ is 14.44-16.80 ng/mL.

Endotoxin

Less than 1.0 EU/ug as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.