



Recombinant Saponaria officinalis Ribosome-inactivating protein saporin-6 (SAP6) & Streptococcus sp. group G Immunoglobulin G-binding protein G (spg)


Product Code	CSB-EP6472
Abbreviation	Recombinant Saponaria officinalis SAP6 & Streptococcus sp spg protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20656&P19909
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Saponaria officinalis Ribosome-inactivating protein saporin-6 (SAP6) & Streptococcus sp. group G Immunoglobulin G-binding protein G (spg)
Immunogen Species	Saponaria officinalis
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	VTSITLDLVNPTAGQYSSFVDKIRNNVKDPNLKYGGTDIAVIGPPSKEKFLRINF QSSRGTVSLGLKRDNLVYVAYLAMDNNTNVNRAYYFRSEITSAESTALFPEATT ANQKALEYTEDYQSIEKNAQITQGDQSRKELGLGIDLLSTSMEAVNKKARVVK DEARFLLIAIQMTAEAAARFRYIQNLVIKNFPNKFENSENKVIQFEVNWKKISTAIYG DAKNGVFNKDYDFGFGKVRQVKDLQMGLLMYLGKPKIDEILAAPKTDITYKLIL NGKTLKGETTTEAVDAATAEKVFKQYANDNGVDGEWYDDATKTFTVTEKPE VIDASELTPAVTTYKLVLINGKTLKGETTTEAVDAATAEKVFKQYANDNGVDGE WYDDATKTFTVTEKPEVIDASELTPAVTTYKLVLINGKTLKGETTTKAVDAETA KAFKQYANDNGVDGVWYDDATKTFTVTE
Research Area	Others
Source	E.coli
Target Names	N/A
Expression Region	25-277aa(P20656?& 291-497aa(P19909)
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged



Mol. Weight 51.8 kDa

Protein Length Linked Heterodimer

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.