




Recombinant Human T-cell differentiation antigen CD6(CD6), partial

Product Code	CSB-YP326953HUa0
Abbreviation	Recombinant Human CD6 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P30203
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Human T-cell differentiation antigen CD6(CD6), partial
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	≥ 95% as determined by SDS-PAGE.
Sequence	HPSPAPPDQLNTSSAESELWEPGERLPVRLTNGSSSCSGTVEVRLEASWEPA CGALWDSRAAEAVCRALGCGGAEAAASQLAPPTPELPPPPAAGNTSVAANATL AGAPALLCSGAEWRLCEVVEHACRS DGRRARVTCAENRALRLVDGGGACAG RVEMLHGEWGSVCDDTDLED AHVVC RQLGCGWAVQALPGLHFTPGRGP IHRDQVNC SGAEAYLWDCPGLPGQHYCGHKEDAGAVCSEHQSWRLTG GAD RCEGQVEVHFRGVWNTVCDSEWYPSEAKVLCQSLGCGTAVERP KGLPHSL S GRMYYS CNGEELT LSNCSWRFNNSNLCSQSLAARVLC SASRSLHNLSTPEVP ASVQTVTIESSVT VKIENKESRE
Research Area	Immunology
Source	Yeast
Target Names	CD6
Expression Region	18-398aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	42.6 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.