





SEZ6 Recombinant Monoclonal Antibody

Product Code	CSB-RA684467MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53EL9
Immunogen	Recombinant Human SEZ6 protein
Species Reactivity	Human, Cynomolgus
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SEZ6
Clone No.	12H9

Image



The Binding Activity of Cynomolgus SEZ6 with Anti-SEZ6 recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Cynomolgus SEZ6 (CSB-MP6419MOV) at 2 µg/mL can bind Anti-SEZ6 recombinant antibody. The EC₅₀ is 0.6865-1.065 ng/mL.



The Binding Activity of Human SEZ6 with Anti-SEZ6 Recombinant Antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human SEZ6 (CSB-MP684467HU2) at 2 µg/mL can bind Anti-SEZ6 recombinant antibody. The EC₅₀ is 2.324-2.770 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.