




GLP1R Recombinant Monoclonal Antibody

Product Code	CSB-RA009514MA3HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43220
Immunogen	Recombinant Human GLP1R protein
Species Reactivity	Mouse
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG4(S228P)
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience;Cardiovascular;Metabolism
Target Names	GLP1R
Clone No.	6C1

Image



The Binding Activity of Mouse Glp1r with Anti-GLP1R Recombinant Antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Mouse Glp1r (CSB-MP009514MO1) at 2 µg/mL can bind Anti-GLP1R recombinant antibody. The EC₅₀ is 141.7-172.4 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.