





PCSK9 Recombinant Monoclonal Antibody

Product Code	CSB-RA017647MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NBP7
Immunogen	Recombinant Human PCSK9 protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer?Cardiovascular?Cell biology;Metabolism;Signal transduction?Stem cells
Target Names	PCSK9
Clone No.	9F1

Image



The Binding Activity of Human PCSK9 with Anti-PCSK9 recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Anti-PCSK9 recombinant antibody at 2 µg/mL can bind Human PCSK9 protein (CSB-MP017647HU(M)-B). The EC₅₀ is 7.450-8.692 ng/mL.



The Binding Activity of Human PCSK9 with Anti-PCSK9 recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human PCSK9 (CSB-MP017647HU(M)-B) at 2 µg/mL on streptavidin coated plates can bind Anti-PCSK9 recombinant antibody. The EC₅₀ is 0.3792-1.157 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.