




Recombinant Human Netrin-1 (NTN1), partial (Active)


Product Code	CSB-MP016127HU(A4)
Abbreviation	Recombinant Human NTN1 protein, partial (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O95631
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm sterile filtered PBS, 6% Trehalose, pH 7.4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Biological Activity	Measured by its binding ability in a functional ELISA. Immobilized Human NTN1 at 2 µg/mL can bind Anti-NTN1 recombinant antibody (CSB-RA016127MA1HU). The EC50 is 26.43-30.77 ng/mL.
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	VRGGPGLSMFAGQAAQPDPDENGHPRRRCIPDFVNAAFGKDVRVSSTCGR PPARYCVVSEERGEERLRSCHLCNASDPKKAHPPAFLDLNNPHNLTCWQSEN YLQFPHNVTLTSLGKKFEVTYVSLQFCSPRPESMAIYKSMDYGRWTWVPFQFY STQCRKMYNRPHRAPITKQNEQEAVCTDSHTDMRPLSGGLIAFSTLDGRPSA HDFDNSPVLQDWVTATDIRVAFSRLHTFGDENEDDSELARDSYFYAVSDLQV GGRCKCNGHAARCVRDRDDSLVCDCRHNTAGPECDRCKPFHYDRPWQRAT AREANECVACNCLHARRCRFMELYKLSGRKSGGVCLNCRHNTAGRHCY CKEGYYRDMGKPITHRKACKACDCHPVGAAGKTCNQTGGQCCKDGVGTGIT CNRCAKGYQQSRSPAPIKIPVAPPTTAASSVEEPEDCDSYCKASKGKLNIM KKYCKKDYAVQIHILKADKAGDWWKFTVNIISVYKQGTSRIRRGDQSLWIRSRD IACKCPKIKPLKYLKLLGNAEDSPDQSGIVADKSSLVIQWRDTRARRLRKFQQR EKKGKCKKA
Source	Mammalian cell
Target Names	NTN1
Expression Region	22-604aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	67.1 kDa
Protein Length	Partial



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Activity
Measured by its binding ability in a functional ELISA. Immobilized Human NTN1 at 2 µg/ml can bind Anti-NTN1 recombinant antibody (CSB-RA016127MA1HU). The EC₅₀ is 26.43-30.77 ng/mL.

Endotoxin

Less than 1.0 EU/ug as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.