





Recombinant *Macaca fascicularis* Natural cytotoxicity triggering receptor 3 (NCR3), partial (Active)

Product Code	CSB-MP015551MOV1
Abbreviation	Recombinant Cynomolgus monkey NCR3 protein, partial (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P61483
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, 6% Trehalose, pH 7.4
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (Cynomolgus monkey)
Biological Activity	Measured by its binding ability in a functional ELISA. Immobilized <i>Macaca fascicularis</i> NCR3 at 2 µg/mL on an Nickel Coated plate can bind <i>Macaca fascicularis</i> NCR3LG1(CSB-MP4942MOVd9). The EC50 is 7.429 to 11.19 ng/mL.
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	LWVSPPEIRTLEGSSAFLPCSFNASQGRLAIGSVTWFRDEVAPGKEVRNGTPEFRGRLAPLSSSRFLRDHQAELHIWDVRGHDAGIYVCRVEVLGGLGVGTGNGTRLVVEKEYPQLG
Source	Mammalian cell
Target Names	NCR3
Expression Region	19-135aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	14.5 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Activity

Measured by its binding ability in a functional ELISA. Immobilized *Macaca fascicularis* NCR3 at 2 µg/ml on an Nickel Coated plate can bind *Macaca fascicularis* NCR3LG1 (CSB-MP4942MOVd9). The EC₅₀ is 7.429 to 11.19 ng/mL.

Endotoxin

Less than 1.0 EU/ug as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.