




Recombinant *Oryza glumipatula* Citrate transporter-like domain-containing protein-VLPs

Product Code	CSB-MP6061GYI
Abbreviation	Recombinant <i>Oryza glumipatula</i> Citrate transporter-like domain-containing protein-VLPs
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0A0D9Z0W3
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, 6% Trehalose, pH 7.4.
Product Type	Recombinant <i>Oryza glumipatula</i> Citrate transporter-like domain-containing protein-VLPs
Immunogen Species	<i>Oryza glumipatula</i>
Sensitivity	Not Test
Purity	VLPs
Sequence	MSELASAPKVALGSIFAVFWMMAVFPSVPFLPIGRTAGSLLSAVLMVIFHVISP DDAYASIDLPLGLLFATMVVGSYLNRNAGMFKHLGRLLAWKSQGGRDLMCRV CVVTALASALFTNDTCCVVLTEFVLELAAERNLPAKPFLALASSANIGSAATPI GNPQNLVIAFNKITFPKFLMGILPAMLVGMVNMVMLLCMYWRELGGGAELS VDGKQMEAVEEGRSPASAKSTPQLNGNGNTMMSLEMSENITTKHPWFMQCT EARRKLFLKSFAYVVTVGMVVAYMVGLNMSWTAITLALVVVDFRDAEPCLD TVSYLLVFFSGMFITVSGFNKTGLPGAIFWDFMAPYSKVN SVGGISVLSVILL SNLASNVPTVLLMGDEVAKAAALISPAAVTTSWLLLAWVSTVAGNLSLLGSAA NLIVCEQARRAPRNAYDLTFWQHIVFGVPSTLIVTAVGIPLIGKI
Research Area	Others
Source	Mammalian cell
Target Names	N/A
Expression Region	1-472aa
Notes	The VLPs are expressed from human 293 cells (HEK293). Mix the sample gently by repeatedly pipetting it up and down. Do not vortex. Repeated freezing and thawing is not recommended. Store the protein at -20°C/-80°C upon receiving it, and ensure to avoid repeated freezing and thawing, otherwise, it will affect the protein activity. The immunization strategy should be optimized (antigen dose, regimen and adjuvant).
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	51.7 kDa

**Protein Length**

Full Length

Image

CSB-MP6061GYI is detected by Mouse anti-6*His monoclonal antibody.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. Aliquot for long-term storage at -80°C. Solubilize for 60 minutes at room temperature with occasional gentle mixing. Avoid vigorous shaking or vortexing.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.