





Recombinant Human Probable G-protein coupled receptor 75 (GPR75)-VLPs

Product Code	CSB-MP009829HU(A4)m8
Abbreviation	Recombinant Human GPR75 protein-VLPs
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O95800
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, 6% Trehalose, pH 7.4.
Product Type	Recombinant Human Probable G-protein coupled receptor 75 (GPR75)-VLPs
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 95% as determined by SEC-HPLC.
Sequence	MNSTGHLQDAPNATSLHVPHSQEGNSTSLQEGQLDIHTATLVCTFLLAVIFC LGSYGNFIVFLSFFDPAFRKFRTNFDLMILNLSFCDLFCGVTAPMFTFVLFSS ASSIPDAFCFTFHLTSSGFIIMSLKTVAVIALHRLRMVLGKQPNRTASFPCTVLLT LLLWATSFTLATLATLKTSKSHLCLPMSSLIAGKKGKAILSLYVVDFTFCVAVVSV SYIMIAQTLRKNAQVRKCPPVITVDASRPQPFMGVPVQGGGDPIQCAMPALYR NQNYNKLQHVQTRGYTKSPNQLVTPAASRLQLVSAINLSTAKDSKAVVTCVIV LSVLVCCPLGISLVQVVLSSNGSFILYQFELFGFTLIFFKSGLNPFYISRNAGL RRKVLWCLQYIGLGFCCCKQKTRLRAMGKGNLEVNRNKSSHETNSAYMLSP KPQKKFVDQACGPSHSKESMVSPKISAGHQHCGQSSSTPINTRIEPYYSIYNS SPSQEESSPCNLQPVNSFGFANSYIAMHYHTTNDLVQEYDSTSAKQIPVPSV
Research Area	Signal Transduction
Source	Mammalian cell
Target Names	GPR75
Expression Region	1-540aa
Notes	The VLPs are expressed from human 293 cells (HEK293). Mix the sample gently by repeatedly pipetting it up and down. Do not vortex. Repeated freezing and thawing is not recommended. Store the protein at -20°C/-80°C upon receiving it, and ensure to avoid repeated freezing and thawing, otherwise, it will affect the protein activity. The immunization strategy should be optimized (antigen dose, regimen and adjuvant).
Tag Info	N-terminal Flag-tagged and C-terminal 6xHis-tagged
Mol. Weight	61.3 kDa


Protein Length
Full Length
Image


CSB-MP009829HU(A4)m8 is detected by Mouse anti-6*His monoclonal antibody.



The purity of VLPs was greater than 95% as determined by SEC-HPLC

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. Aliquot for long-term storage at -80°C. Solubilize for 60 minutes at room temperature with occasional gentle mixing. Avoid vigorous shaking or vortexing.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.