





Recombinant Dog Programmed Cell Death Protein 1 precursor, partial, Biotinylated

Product Code	CSB-MP5072DOj8-B
Abbreviation	Recombinant Dog Programmed Cell Death protein 1 precursor protein, partial, Biotinylated
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	NP_001301026
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Canis familiaris Programmed Cell Death Protein 1 precursor, partial
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Sensitivity	Not Test
Purity	≥ 90% as determined by SDS-PAGE. ≥ 95% as determined by SEC-HPLC.
Sequence	LDSPDRPWSPLTFSPAQLTVQEGENATFTCSLADIPDSFVLNWyRLSPRNQTD KLAAFQEDRIEPGRDRRFRVTRLPNGRDFHMSIVAARLNDSGIYLCGAIYLPPN TQINESPRAELSVTERTLEPPTQSPSPPPRLSGQLQG
Research Area	Others
Source	Mammalian cell
Target Names	N/A
Expression Region	25-168aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal hFc1-Avi-tagged
Mol. Weight	44.0 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



The purity of 5072 was greater than 95% as determined by SEC-HPLC

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.