



# Recombinant *Macaca fascicularis* Heat shock protein HSP 90-beta (HSP90AB1)

<b>Product Code</b>	CSB-MP686992MOV
<b>Abbreviation</b>	Recombinant <i>Cynomolgus</i> monkey HSP90AB1 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q4R4T5
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
<b>Product Type</b>	Recombinant <i>Macaca fascicularis</i> Heat shock protein HSP 90-beta (HSP90AB1)
<b>Immunogen Species</b>	<i>Macaca fascicularis</i> (Crab-eating macaque) ( <i>Cynomolgus</i> monkey)
<b>Sensitivity</b>	Not Test
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	MPEEVHHGEEEVETFAFQAEIAQLMSLIINTFYSNKEIFLRELISNASDALDKIRY ESLTDPSKLD SGKELKIDIIPNPQERTLTLVDTGIGMTKADLNNLGTIAKSGTKA FMEALQAGADISMIGQFGVGFYSAYLVAEKVVVITKHNDDEQYAWESSAGGSF TVRADHGEPIGRGTKVILHLKEDQTEYLEERRVKEVVKKHSQFIGYPITLYLEKE REKEISDDEAEKKEEEDKDDEEKPKIEDVGSDEEDDSGKDKKKKTKKIK EKYIDQEELNKTPIWTRNPDDITQEEYGEFYKSLTNDWEDHLAVKHFSVEGQ LEFRALLFIPRRAPFDLFENKKKNNIKLYVRRVFIMDSCDELIPEYLNFI RGVVD SEDLPLNISREMLQQSKILKVIRKNIVKKCLELFS ELAEDKENYKKFYEA FSKNL KLG IHEDSTNRRRLSELLRYHTSQSGDEMTSLSEYVSRMKETQKSIYYITGESK EQVANS AFVERVRKRGFEVVMTEPID EYCVQQLKEFDGKSLVSVTKEGLELP EDEEEKKMEESKAKFENLCKLMKEILDKKVEKVTISNRLVSSPCCIVTSTYGW TANMERIMKAQALRDNSTMGYMMAKKHLEINPDHPIVETLRQKAEADKNDKAV KDLVLLFETALLSSGFSLEDPQTHSNRIYRMIKLG LGIDEDEVA AEPSAAVPD EIPPLEGDE D ASRMEEVD
<b>Research Area</b>	Others
<b>Source</b>	Mammalian cell
<b>Target Names</b>	HSP90AB1
<b>Expression Region</b>	1-724aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	C-terminal 10xHis-tagged



---


**Mol. Weight** 84.7 kDa

---

**Protein Length** Full Length

---

**Image**



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.