



NLRP12 Antibody

Product Code	CSB-PA015867GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P59046
Immunogen	Recombinant Human NLRP12
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IP
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NLRP12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.