




Recombinant Plutella xylostella N6-adenosine-methyltransferase subunit METTL14

Product Code	CSB-BP6489EXE
Abbreviation	Recombinant Plutella xylostella N6-adenosine-methyltransferase subunit METTL14 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	WQH20074.1
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Plutella xylostella N6-adenosine-methyltransferase subunit METTL14
Immunogen Species	Plutella xylostella (Diamondback moth) (Plutella maculipennis)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MSEKLKELRERSQKRKKLLAQT LGVSSVSEL RHAL GASLDVPAKKQAAAEGG PDDDKPAAKDQPD ELVYTD SSTFLKGTQSSNPHNDYCQHFVDTGQRPQNFIR DVGLADRFE EYPKLREL IKLKDD LIAETATPPMYLKCDLKT FDLKELNNKF DVLL VEPPLGGAWRWRD VLDLELQQLAQPRAFVFLWCGSSDGLDMGRECLKK WG FRRCEDICWIKTNIN NPGH SKNLEQNAV FQRTKEHCLMG IKGTVRRSVDGDFI HANVDIDL ISEDAEFGSAEKPIE IFHIMEHFCLGRRRIHLFGRDSTIRPGWVTLG HELTNSNFNASLYASYFLGGRSHTGCTERIEALRPKSPPAAGKTRPRRPGFRP RPRTNV
Research Area	Others
Source	Baculovirus
Target Names	N/A
Expression Region	1-377aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	43.6 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.