



Recombinant Human Alpha-fetoprotein (AFP)

| | |
|--------------------------|--|
| Product Code | CSB-EP001421HU |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P02771 |
| Form | Liquid or Lyophilized powder |
| Storage Buffer | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0. |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | RTLHRNEYGIASILDSYQCTAEISLADLATIFFAQFVQEATYKEVSKMVKDALTAI EKPTGDEQSSGCLLENQLPAFLEELCHEKEILEKYGHSDCCSQSEEGRHNCFL AHKKPTPASIPLFQVPEPVTSCAEYEEDRETFMKNFIYEIARRHPFLYAPTILLW AARYDKIIPSCCKAENAVECFQTKAATVTKELRESSLLNQHACAVMKNFGTRT FQAITVTKLSQKFTKVNFTIEIQKLVLDVAHVHEHCCRGDVLDCLDQDGEKIMSYI CSQQDTLSNKITECCKLTTLARGQCIIHAENDEKPEGLSPNLNRFLGDRDFNQF SSGEKNIFLASFVHEYSRRHPQLAVSVILRVAKGYQELLEKCFQTENPLECQDK GEEELQKYIQESQALAKRSCGLFQKLGEYYLQNAFLVAYTKKAPQLTSSELMAI TRKMAATAATCCQLSEDKLLACGEGAADIIIGHLCIRHEMTPVNPVGVGQCCTSS YANRRPCFSSLVVDETYVPPAFSDDKFIFHKDLCQAQGVALQTMKQEFNLNV KQKPQITEEQLEAVIADFSGLLEKCCQGQEQEVCFEAEEGQKLIKTRAALGV |
| Source | E.coli |
| Target Names | AFP |
| Expression Region | 19-609aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.