




CYC1 Recombinant Monoclonal Antibody

Product Code	CSB-RA006327MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08574
Immunogen	Recombinant Human CYC1 protein
Species Reactivity	Human
Tested Applications	ELISA, WB, FC; Recommended dilution: WB:1:500-1:2000, FC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYC1
Clone No.	27E12D11

Image



Western Blot

Positive WB detected in: A549 whole cell lysate(20μg), PC-3 whole cell lysate(20μg), HeLa whole cell lysate(20μg), K562 whole cell lysate(20μg), HepG2 whole cell lysate(20μg), Jurkat whole cell lysate(20μg), SH-SY5Y whole cell lysate(20μg)

All lanes: CYC1 antibody at 1:1000


Secondary

Goat polyclonal to human IgG at 1/40000 dilution

Predicted band size: 35 kDa

Observed band size: 35 kDa

Exposure time?1min



Overlay Peak curve showing PC-3 cells stained with CSB-RA006327MA1HU (red line) at 1:100. The cells were fixed in 4% formaldehyde and permeated by 0.2% TritonX-100 for 10min. Then 10% normal goat serum to block non-specific protein-protein interactions followed by the antibody (1 μ g/1*10⁶cells) for 45min at 4?. The secondary antibody used was Fluorescein (FITC) AffiniPure Goat Anti-Human IgG, Fc γ fragment specific at 1:200 dilution for 35min at 4?. Control antibody (green line) was human IgG1 (1 μ g/1*10⁶cells) used under the same conditions. Acquisition of >10,000 events was performed.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.