




# PDCD1LG2 Recombinant Monoclonal Antibody

<b>Product Code</b>	CSB-RA017667MA1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9BQ51
<b>Immunogen</b>	Recombinant Human IL17A protein
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	Affinity-chromatography
<b>Isotype</b>	hIgG1
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Recombinant Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	IL17A Recombinant Monoclonal Antibody
<b>Clone No.</b>	7F11D11

## Image



### Western Blot

Positive WB detected in: HepG2 whole cell lysate(20µg), HeLa whole cell lysate(20µg), A549 whole cell lysate(20µg)

All lanes: PDCD1LG2 antibody at 1:1000

Secondary

Goat polyclonal to human IgG at 1/40000 dilution

Predicted band size: 31 kDa

Observed band size: 50 kDa

Exposure time?2min

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.