





# SERPINE1 Recombinant Monoclonal Antibody

<b>Product Code</b>	CSB-RA021081MA1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P05121
<b>Immunogen</b>	Recombinant Human SERPINE1 protein
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, FC; Recommended dilution: WB:1:500-1:2000, FC:1:50-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	Affinity-chromatography
<b>Isotype</b>	hIgG1
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Recombinant Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	SERPINE1
<b>Clone No.</b>	9H8

## Image



Western Blot  
 Positive WB detected in: U-251MG whole cell lysate(20µg)  
 All lanes: SERPINE1 antibody at 1:1000  
 Secondary  
 Goat polyclonal to human IgG at 1/40000 dilution  
 Predicted band size: 45 kDa  
 Observed band size: 45 kDa  
 Exposure time?5min



Overlay Peak curve showing 786-O cells stained with CSB-RA021081MA1HU (red line) at 1:100. The cells were fixed in 4% formaldehyde and permeated by 0.2% TritonX-100 for 10min. Then 10% normal goat serum to block non-specific protein-protein interactions followed by the antibody (1ug/1\*10<sup>6</sup>cells) for 45min at 4?. The secondary antibody used was Fluorescein (FITC) AffiniPure Goat Anti-Human IgG, Fc $\gamma$  fragment specific at 1:200 dilution for 35min at 4?. Control antibody (green line) was human IgG1 (1ug/1\*10<sup>6</sup>cells) used under the same conditions. Acquisition of >10,000 events was performed.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.