




FCER1A Recombinant Monoclonal Antibody

Product Code	CSB-RA008532MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12319
Immunogen	Recombinant Human FCER1A protein
Species Reactivity	Human
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:500-1:2000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FCER1A
Clone No.	15E5

Image



Western Blot

Positive WB detected in: K562 whole cell lysate(20μg)

All lanes: FCER1A antibody at 1:1000


Secondary

Goat polyclonal to human IgG at 1/40000 dilution

Predicted band size: 30 kDa

Observed band size: 30 kDa

Exposure time?5min



Immunofluorescence staining of HeLa cell with CSB-RA008532MA1HU at 1:30 ,counter-stained with DAPI. The cells were fixed in 4% formaldehyde and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Fluorescein (FITC) AffiniPure Goat Anti-Human IgG, Fcγ fragment specific.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.