




NAT10 Recombinant Monoclonal Antibody

Product Code	CSB-RA189183A0HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0A0
Immunogen	A synthesized peptide derived from human NAT10
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, FC; Recommended dilution: WB:1:500-1:2000, FC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Purification Method	Affinity-chromatography
Isotype	Rabbit IgG
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NAT10
Clone No.	9C11

Image




Western Blot

Positive WB detected in: SH-SY5Y whole cell lysate(20µg), A549 whole cell lysate(20µg), HeLa whole cell lysate(20µg), NIH/3T3 whole cell lysate(20µg), HEK293 whole cell lysate(20µg), PC-3 whole cell lysate(20µg), COLO205 whole cell lysate(20µg), Jurkat whole cell lysate(20µg), HepG2 whole cell lysate(20µg)
All lanes: NAT10 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/40000 dilution
Predicted band size: 116 kDa
Observed band size: 130 kDa
Exposure time?20s



Overlay Peak curve showing jurkat cells stained with CSB-RA189183A0HU (red line) at 1:100. The cells were fixed in 4% formaldehyde and permeated by 0.2% TritonX-100 for 10min. Then 10% normal goat serum to block non-specific protein-protein interactions followed by the antibody (1ug/1*10⁶cells) for 45min at 4?. The secondary antibody used was FITC-conjugated goat anti-rabbit IgG (H+L) at 1/200 dilution for 35min at 4?.Control antibody (green line) was Rabbit IgG (1ug/1*10⁶cells) used under the same conditions. Acquisition of >10,000 events was performed.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.