




TUBB4B/TUBB2A/TUBB2B Recombinant Monoclonal Antibody

Product Code	CSB-RA130268A0HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68371/Q13885/Q9BVA1
Immunogen	A synthesized peptide derived from human TUBB4B/TUBB2A/TUBB2B
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, FC; Recommended dilution: WB:1:500-1:2000, FC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Purification Method	Affinity-chromatography
Isotype	Rabbit IgG
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TUBB4B/TUBB2A/TUBB2B
Clone No.	5A9

Image



Western Blot

Positive WB detected in: SH-SY5Y whole cell lysate(30µg), Mouse brain tissue lysate(30µg), Rat brain tissue lysate(30µg)

All lanes: beta 4B/2A/2B Tubulin antibody at 1:1000


Secondary

Goat polyclonal to rabbit IgG at 1/40000 dilution

Predicted band size: 50 kDa

Observed band size: 55 kDa

Exposure time?40s



Overlay Peak curve showing HeLa cells stained with CSB-RA130268A0HU (red line) at 1:100. The cells were fixed in 4% formaldehyde and permeated by 0.2% TritonX-100 for 10min. Then 10% normal goat serum to block non-specific protein-protein interactions followed by the antibody (1 μ g/1 $\times$ 10⁶cells) for 45min at 4 $^{\circ}$ C. The secondary antibody used was FITC-conjugated goat anti-rabbit IgG (H+L) at 1/200 dilution for 35min at 4 $^{\circ}$ C. Control antibody (green line) was Rabbit IgG (1 μ g/1 $\times$ 10⁶cells) used under the same conditions. Acquisition of >10,000 events was performed.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.