



# Recombinant Mouse Sarcoplasmic/endoplasmic reticulum calcium ATPase 3 (Atp2a3), partial

<b>Product Code</b>	CSB-EP727334MO
<b>Abbreviation</b>	Recombinant Mouse Atp2a3 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q64518
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
<b>Product Type</b>	Recombinant Mouse Sarcoplasmic/endoplasmic reticulum calcium ATPase 3 (Atp2a3), partial
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Sensitivity</b>	Not Test
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	VITTCLALGTRRMARKNAIVRSLPSVETLGCTSVICSDKTGTLTTNQMSVCRMF VVAEAEAGTCRLHEFTISGTTYTPEGEVRQGEQPVRCGQFDGLVELATICALC NDSALDYNEAKGVYEKVGAEATETALTCLVEKMNVFDTDLKGLSRVERAGACN SVIKQLMRKEFTLEFSRDRKSMVYCTPTRADPKVQGSKMFVKGAPESVIERC SSVRVGSRTAPLSTTSREHILAKIRDWGS S DTLRCLALATRDTPPRKEDMHL DDCSRFBVQYETDLTFVGCVGM L D P P R P E V A A C I T R C S R A G I R V V M I T G D N K G T A V A I C R R L G I F G D T E D V L G K A Y T G R E F D D L S P E Q Q R Q A C R T A R C F A R V E P A H K S R I V E N L Q S F N E I T A M T G D G V N D A P A L K K A E I G I A M G S G T A V A K S A A E M V L S D D N F A S I V A A V E E G R A I Y N N M
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	Atp2a3
<b>Expression Region</b>	314-757aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	C-terminal 6xHis-tagged
<b>Mol. Weight</b>	55.1 kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.