




Recombinant Pongo abelii Dyslexia-associated protein KIAA0319-like protein (KIAA0319L), partial

Product Code	CSB-YP735876PYX
Abbreviation	Recombinant Pongo abelii Kiaa0319l protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RFR6
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Pongo abelii Dyslexia-associated protein KIAA0319-like protein (KIAA0319L), partial
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	VVSAGDSVQITLPKNEVQLNAYVLQEPPKGETYTYDWQLITHPRDYSGEMEG KHSQILKLSKLTPLGYEFKVIVEGQNAHGEGYVNVTVKPEPRKNRPPIAIVSPQ FQEISLPTTSTVIDGSQSTDDDKIVQYHWHEELKGPLREEKISEDTAILKLSKLVP GNYTFSLTVVSDGATNSTTANLTVNKAVDYPPVANAGPNQVITLPQNSITLFG NQSTDDHGITSYEWSLSPSSKGKVVEMQGVRTPTLQLSAMQEGDYTYQLTVT DTIGQQATAQVTVIVQPENNKPPQADAGPDKELTLPVDSTTLTGSKSSDDQKII SYLWEKTQGPDGVQLENANSSIATVTGLQVGTYVFTLTVKDERNLQSQSSVN VIVKEEINKPPIAKITGNVVITLPTSTAELDGSKSSDDKGIVSYLWTRDEGSPAA GEVLNHSDDHHPILFLSNLVEGTYTFHLKVTDAGGESDTRTTVEVKPD
Research Area	Others
Source	Yeast
Target Names	Kiaa0319l
Expression Region	310-785aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	53.3 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.