




# Recombinant Strongyloides stercoralis L3NieAg.01 (L3Nie.01) (X208Y,X221W), partial

<b>Product Code</b>	CSB-EP4964GLF
<b>Abbreviation</b>	Recombinant Strongyloides stercoralis L3Nie.01 protein (X208Y,X221W), partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UA16
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
<b>Product Type</b>	Recombinant Strongyloides stercoralis L3NieAg.01 (L3Nie.01) (X208Y,X221W), partial
<b>Immunogen Species</b>	Strongyloides stercoralis (Threadworm)
<b>Sensitivity</b>	Not Test
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	NSARDENGKLIYTYNGNDYDTKEAMEDAIQRDYPDKIFTFGGDNNNGKKRKIDI SKWKGNNTFSNKIFDEIWEGYNYDNDKAKNFKVMKTKLFNEQNKYRIAHGAK KLIKSKDLEKKAQAYAEVIARLGRLEHDPKNRIEGTGENLAYGTTFIGHLAVKG WYDEIALYNFKKPGFSPATGHFTQLVWKGTTTHAGFGVVEKGDREVYVYKYSP PGNYPRQFWANVLQRKQ
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	L3Nie.01
<b>Expression Region</b>	1-229aa(X208Y,X221W)
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	C-terminal 6xHis-tagged
<b>Mol. Weight</b>	33.3 kDa
<b>Protein Length</b>	Partial
<b>Image</b>	 <p>(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.</p>

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.