




SLC3A2 Recombinant Monoclonal Antibody

Product Code	CSB-RA021640MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08195
Immunogen	Recombinant Human SLC3A2 protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC3A2
Clone No.	2B8

Image



The Binding Activity of Human SLC3A2 with Anti-SLC3A2 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human SLC3A2 protein (CSB-MP021640HU1) at 2 µg/mL can bind Anti-SLC3A2 recombinant antibody. The EC₅₀ is 2.257-3.091 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.