




Recombinant Staphylococcus epidermidis Ornithine carbamoyltransferase (argF)

Product Code	CSB-EP017270SLD
Abbreviation	Recombinant Staphylococcus epidermidis argF protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C0N1
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Staphylococcus epidermidis Ornithine carbamoyltransferase (argF)
Immunogen Species	Staphylococcus epidermidis (strain ATCC 12228 / FDA PCI 1200)
Sensitivity	Not Test
Purity	≥ 95% as determined by SDS-PAGE.
Sequence	MKNLRNRSFLTLLDFSRQEVEFLLTSEDLKRAKYIGTEKPMMLKNKNIALLFKED STRTRCAFEVAAHQGAHVTYLGPTGSQMGGKETAKDTARVLGGMYDGIEYR GFSQRTVETLAQYSGVPVWNGLTDEDHPTQVLADFLTAKVLLKKEYADINFTY VGDRNNVANALMQGAAIMGMNFHLVCPKELNPTEELLNRCERIATENGGNIL ITDDIDKGVKDSVDIYTDVWVSMGEPDEVWQERLKLKPYQVNQALLEKTGNP NVIFEHCLPSFHNAETKIGQIQIYEKYGISEMEVTDDVFESKASVVFQEAENRMH TIKAVMVATLGEF
Research Area	Others
Source	E.coli
Target Names	argF
Expression Region	1-333aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	44.5 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.