




Recombinant Human Biotin--protein ligase (HLCS)

Product Code	CSB-EP010518HU
Abbreviation	Recombinant Human HLCS protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P50747
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Human Biotin--protein ligase (HLCS)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than $\geq 85\%$ as determined by SDS-PAGE.
Sequence	MEDRLHMDNGLVPQKIVSVHLQDSTLKEVKDQVSNKQAQILEPKPEPSLEIKP EQDGMHEHVGRDDPKALGEEPQRRGSASGSEPAGDSDRGGGPVEHYHLHL SSCHECLELENSTIESVKFASAENIPDLPYDYSSSLESVADETSPEREGRRVNL TGKAPNILLYVGSQSQEALGRFHEVRSVLADCVDIDSYILYHLLSDSALRDPWT DNCLLLVIATRESIPEDLYQKFMAYLSQGGKVLGLSSSFTFGGFQVTSKGALHK TVQNLVFSKADQSEVKLSVLSSGCRYQEGPVRLSPGRLQGHLNEDKDRMIV HVPFGTRGGEAVLCQVHLELPPSSNIVQTPEDFNLLKSSNFRRYEVLREILTTL GLSCDMKQVPALTPLYLLSAAEEIRDPLMQWLGKHVDSEGEIKSGQLSLRFVS SYVSEVEITPSCIPVVTNMEAFSSEHFNLEIYRQNLQTKQLGKVVILFAEVTPTTM RLLDGLMFQTPQEMGLIVIAARQTEGKGRGGNVWLSPVGCALSTLLISIPRSQ LGQRIPFVQHLMSVAVVEAVRSIPEYQDINLRVKWPNDIYSDLMKIGGVLVNS TLMGETFYILIGCGFNVTNSNPTICINDLITEYNKQHKAEKPLRADYLIARVVTV LEKLIKEFQDKGPNVSLPLYRYWVHSGQQVHLGSAEGPKVSIVGLDDSGFLQ VHQEGGEVTVHPDGNSFDMLRNLIIPKRR
Research Area	Signal Transduction
Source	E.coli
Target Names	HLCS
Expression Region	1-726aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	85.8 kDa

**Protein Length****Full Length****Image**

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.