



Recombinant Rat Heat shock protein HSP 90-alpha (Hsp90aa1)


Product Code	CSB-EP010802RA
Abbreviation	Recombinant Rat Hsp90aa1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P82995
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Rat Heat shock protein HSP 90-alpha (Hsp90aa1)
Immunogen Species	Rattus norvegicus (Rat)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	<p>PEETQTQDQPMEEEEVETFAFQAEIAQLMSLIINTFYSNKEIFLRELISNSSDAL DKIRYESLTDPSKLDGKELHINLIPNKQDRTLTIVDTGIGMTKADLNNLGTIAK SGTKAFMEALQAGADISMIGQFGVGFYSAYLVAEKVTVITKHNDDEQYAWESS AGGSFTVRTDTGPEPMGRGTKVILHLKEDQTEYLEERRIKEIVKKHSQFIGYPITL FVEKERDKEVSDDEAEKEEKEEKEEKEEKEEESDDKPEIEDVGSDEEEEEKKD GDKKKKKIKIKEYIDQEELNKTPIWTRNPDDITNEEYGEFYKSLTNDWEEHLA VKHFSVEGQLEFRALLFVPRRAPFDLFENRKKKNNIKLYVRRVFIMDNCEELIP EYLNfirGVVdSEDLPLNISREMLQQSKILKVirKnlVkkCLELftelaEDKENY KKfYEQfSKNIKlGIHEDSQNRKKLSELLRYyTSASGDEMvSLKDYCTRmKEN QKHIFITGETKDQVANSaFVERLRKHGLEVIYMIePIDEYCVQQLKEfEGKTLV SVTKEGLELPEDEEEKKKQEEKTKfENLCKIMKDILEKKVekVvVSNRLVTSP CCIVTSTYGWTANMERIMKAQALRDNSTMGYMAAKKHLEINPDHSIIETLRQKA EADKNDKSVKDLVILLYETALLSSGFSLEDPQTHANRIYRMIKLGLGIDEDDPTV DDTSAAVTEEMPpLEGGDDDTSRMEEVD</p>
Research Area	Signal Transduction
Source	E.coli
Target Names	Hsp90aa1
Expression Region	2-733aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged

**Mol. Weight**

91.6 kDa

Protein Length

Full Length of Mature Protein

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.