




Recombinant Triticum aestivum Imidazoleglycerol-phosphate dehydratase 1, chloroplastic

Product Code	CSB-EP335890TQN
Abbreviation	Recombinant Triticum aestivum CAMPLR22A2D_LOCUS4590 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34048
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Triticum aestivum Imidazoleglycerol-phosphate dehydratase 1, chloroplastic
Immunogen Species	Triticum aestivum (Wheat)
Sensitivity	Not Test
Purity	≥ 95% as determined by SDS-PAGE.
Sequence	ASGVGGNGSPMAPEESTVSSRLGEVKRVTKETNVHVKINLDGTGVANSSTGIP FLDHMLDQLASHGLFDVYVKATGDTHIDHHSNEDIALAIGTALLQALGDRKGI NRFGHFTAPLDEAAVEVILDLSGRPHLSCGLSIPTERVGTYDTQLVEHFFQSLV NTSGMTLHIRQLAGNNSHHIEATFKAFARALRQATEYDLRRQGTIPSSKGVLS RS
Research Area	Others
Source	E.coli
Target Names	N/A
Expression Region	53-269aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	30.2 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.