



# CCR4 Monoclonal Antibody

<b>Product Code</b>	CSB-MA236376
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P51679
<b>Immunogen</b>	Human CCR4 transfected cells
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Homo sapiens (Human)
<b>Tested Applications</b>	ELISA, FC; Recommended dilution: FC:1:50-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide
<b>Purification Method</b>	The antibody was purified by affinity chromatography.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Monoclonal Antibody
<b>Target Names</b>	CCR4
<b>Clone No.</b>	1H4
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.