



# 6\*His Monoclonal Antibody, Biotin conjugated

<b>Product Code</b>	CSB-MA000011D0m
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	The immunogen is a mixture of various His-tagged proteins
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	All
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Tag Control Antibody
<b>Clone No.</b>	3G52H8
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.