





FAP Recombinant Monoclonal Antibody

Product Code	CSB-RA008424MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12884
Immunogen	Recombinant Human FAP protein
Species Reactivity	Human, Cynomolgus
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FAP
Clone No.	9B5

Image



The Binding Activity of Cynomolgus FAP with Anti-FAP Recombinant Antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Cynomolgus FAP (CSB-MP6769MOV) at 2 µg/mL can bind Anti-FAP recombinant antibody. The EC₅₀ is 0.6459-0.8172ng/mL.



The Binding Activity of Human FAP with Anti-FAP Recombinant Antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human FAP (CSB-MP008424HU) at 2 µg/mL can bind Anti-FAP recombinant antibody. The EC₅₀ is 0.7317-0.9126 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.