




SDCBP Recombinant Monoclonal Antibody

Product Code	CSB-RA020893MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00560
Immunogen	Recombinant Human SDCBP protein
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SDCBP
Clone No.	8C9

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate(20µg), A549 whole cell lysate(20µg), SH-SY5Y whole cell lysate(20µg), colo205 whole cell lysate(20µg), HepG2 whole cell lysate(20µg), PC-3 whole cell lysate(20µg)

All lanes: SDCBP antibody at 1:1000

Secondary

Goat polyclonal to human IgG at 1/40000 dilution

Predicted band size: 32 kDa

Observed band size: 32 kDa

Exposure time: 2min

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.