




TAF15 Recombinant Monoclonal Antibody

Product Code	CSB-RA899164A0HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92804
Immunogen	A synthesized peptide from human TAF15 protein
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF, FC; Recommended dilution: WB:1:500-1:2000, IHC:1:50-1:200, IF:1:50-1:200, FC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Purification Method	Affinity-chromatography
Isotype	Rabbit IgG
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TAF15
Clone No.	4F11

Image



Western Blot

Positive WB detected in:A431 HeLa whole cell lysate(30µg)

All lanes: TAF15 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/40000 dilution
10,000 events was performed.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.