



Dper\GL27344 Antibody

Product Code	CSB-PA6488XA01DMB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B4GZD1
Immunogen	Synthesized peptide derived from internal of Drosophila persimilis (Fruit fly) DperGL27344 protein (182-201aa)
Raised In	Rabbit
Species Reactivity	Drosophila persimilis (Fruit fly)
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila persimilis (Fruit fly)
Target Names	DperGL27344
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.