



Phospho-CENPC (S261) Antibody

Product Code	CSB-PA005207XA11phHU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03188
Immunogen	Synthesized peptide derived from internal of Homo sapiens (Human) CENPC protein (252-266aa)
Raised In	Rabbit
Species Reactivity	Homo sapiens (Human)
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CENPC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.