



# Recombinant Human Phospholipase A-2-activating protein (PLAA)

<b>Product Code</b>	CSB-YP018107HU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y263
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTSGATRYRLSCSLRGHELDVRGLVCCAYPPGAFVSVSRDRTTRLWAPDSPN RSFTEMHCMMSGHSNFVSCVCIIPSSDIYPHGLIATGGNDHNICIFSLDSPMPLYI LKGHKNTVCSLSSGKFGTLLSGSWDTTAKVWLNDKCMMTLQGHTAAVWAVKI LPEQGLMLTGSADKTVKLWKAGRCERTFSGHEDCVRGLAILSETFLSCANDA SIRRWQITGECLEVYYGHTNYIYSISVFPNCRDFVTTAEDRSLRIWKHGECAQT IRLPAQSIWCCCVDNGDIVVGASDGIIRVFTESEDRTASAEIKAFAKELSHATI DSKTGDLGDINAEQLPGREHLNEPGTREGQTRLIRDGEKVEAYQWSVSEGR WIKIGDVVGSSGANQQTSGKVLVEGKEFDYVFSIDVNEGGPSYKLPYNTSDDP WLTAYNFLQKNDLNPMFLDQVAKFIIDNTKGQMLGLGNPSFSDPFTGGGRYV PGSSGSSNTLPTADPFTGAGRYVPGSASMGTTMAGVDPFTGNSAYRSAASK TMNIYFPKKEAVTFDQANPTQILGKLKELNGTAPEEKKLTEDDLILLEKILSLICN SSSEKPTVQQQLILWKAINCPEDIVFPALDILRLSIKHPSVNENFCNEKEGAQFS SHLINLLNPKGKPANQLLALRTFCNCFVGQAGQKLMMSQRESLMSHAIELKSG SNKNIHIALATLALNYSVCFHKDHNIEGKAQCLSLISTILEVVQDLEATFRLLVAL GTLISDDSNVQLAKSLGVDSQIKKYSSVSEPAKVSECCRFILNLL
<b>Source</b>	Yeast
<b>Target Names</b>	PLAA
<b>Expression Region</b>	1-795
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.