




Recombinant Human C-X-C chemokine receptor type 2 (CXCR2)-VLPs (Active)


Product Code	CSB-MP011673HU(A4)
Abbreviation	Recombinant Human CXCR2 protein-VLPs (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25025
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm sterile filtered PBS, 6% Trehalose, pH 7.4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Biological Activity	Measured by its binding ability in a functional ELISA. Immobilized Human CXCR2 at 2 µg/mL can bind Anti-CXCR2 recombinant antibody (CSB-RA011673MA1HU). The EC50 is 0.3741-0.4671 ng/mL. The VLPs (CSB-MP3838) is negative control.
Purity	VLPs
Sequence	MEDFNMESDSFEDFWKGEDLSNYSYSSTLPPFLLDAAPCEPESLEINKYFVVII YALVFLLSLLGNSLVMVLVILYSRVGRSVTDVYLLNLALADLLFALTLPIWAASKV NGWIFGTFLCKVVSLLEKVNIFYSGILLACISVDRYLAIVHATRRLTQKRYLVKFI CLSIWGLSLLLALPVLLFRRTVYSSNVSPACYEDMGNNNTANWRMLLRILPQSF GFIVPLLIMLFCYGFTRLRFLKAHMGQKHRAMRVIFAVVLIFLLCWLPYNLVLLA DTLMRTQVIQETCERRNHIDRALDATEILGILHSCLNPLIYAFIQKFRHGLLKIL AIHGLISKDSLKPSRPSFVGSSTGHTSTTL
Source	Mammalian cell
Target Names	CXCR2
Expression Region	1-360aa
Notes	The VLPs are expressed from human 293 cells (HEK293). Mix the sample gently by repeatedly pipetting it up and down. Do not vortex. Repeated freezing and thawing is not recommended. Store the protein at -20°C/-80°C upon receiving it, and ensure to avoid repeated freezing and thawing, otherwise, it will affect the protein activity. The immunization strategy should be optimized (antigen dose, regimen and adjuvant).
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	42.2 kDa
Protein Length	Full Length



Image



CSB-MP011673HU(A4) is detected by Mouse anti-6*His monoclonal antibody.



Activity

Measured by its binding ability in a functional ELISA. Immobilized Human CXCR2 at 2 µg/ml can bind Anti-CXCR2 recombinant antibody (CSB-RA011673MA1HU). The EC₅₀ is 0.3741-0.4671 ng/mL. The VLPs (CSB-MP3838) is negative control.

Endotoxin

Less than 1.0 EU/ug as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.