




Recombinant Human Large neutral amino acids transporter small subunit 2 (SLC7A8)

Product Code	CSB-CF887037HU(A4)
Abbreviation	Recombinant Human SLC7A8 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UHI5
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Human Large neutral amino acids transporter small subunit 2 (SLC7A8)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MEEGARHRNNTTEKKHPGGGESDASPEAGSGGGGVALKKEIGLVSACGIIVGNI IGSGIFVSPKGVLENAGSVGLALIVWIVTGFITVVGALCYAELGVTIPKSGGDYS YVKDIFGGLAGFLRLWIAVLVIYPTNQAVIALTFSNYVLQPLFPTCFPPESGLRLL AAICLLLLTWVNCSSVRWATRVQDIFTAGKLLALALIIIMGIVQICKGEYFWLEPK NAFENFQEPDIGLVALAFLQGSFAYGGWNFLNYVTEELVDPYKNLPRAIFISIPL VTFVYVFANVAYVTAMSPQELLASNAVAVTFGEKLLGVMAMWIMPISVALSTFG GVNGSLFTSSRLFFAGAREGHLPSVLAMIHVKRCTPIPALLFTCISTLLMLVTSD MYTLINYVGFINYLFYGVTVAGQIVLRWKKPDIPRIKINLLFPIIYLLFWAFLLVF SLWSEPVVCGIGLAIMLTGVPVYFLGVYWQHKKPKCFSDFIELLTLVSQKMCVV VYPEVERGSGTEEANEDMEEQQQPMYQPTPTKDKDVAGQPQP
Research Area	Neuroscience
Source	in vitro E.coli expression system
Target Names	SLC7A8
Expression Region	1-535aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Twin-Strep-tagged
Mol. Weight	62.9 kDa
Protein Length	Full Length



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.