



Recombinant Human Dual specificity calcium/calmodulin-dependent 3',5'-cyclic nucleotide phosphodiesterase 1A (PDE1A)

Product Code	CSB-EP2273HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P54750
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGSSATEIEELENTTFKYL TGEQTEKMWQRLKGILRCLVKQLERGDVNVVDLK KNIEYAASVLEAVYIDETRLLDTEDELSDIQTDSVPSEVRDWLASTFTRKMG TKKKPEEKPKFRSIVHAVQAGIFVERMYRKTYH MVGLAYPAAVIVTLKDVDKW SFDVFALNEASGEHSLKFM IYELFTRYDLINRFKIPVSCLITFAEALVGY SKYKN PYHNL IHAADV TQTVHYIMLHTGIMHWL TELEILAMVFAAAIHDYEHTGTTN NFH IQTRSDVAILYND RSVLENHHVSAAYRLMQEEEMNILINLSKDDWRDLRNLVIE MVLSTDM SGHFQQIKNIRNSLQQPEGIDRAKTMSLILHAADISHPAKSWKLHYR WTMALMEEFFLQGDKEAELGLPF SPLCDRKSTMVAQSQIGFIDFIVEPTF SLLT DSTEKIVIPLIEEASKAETSSYVASSSTTIVGLHIADALRRSNTKGSMSDGSYSP DYSLAAVDLKSFKNNLVDIIQQNKERWKELAAQEARTSSQKCEFIHQ
Source	E.coli
Target Names	PDE1A
Expression Region	1-535
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.