



# Recombinant Human Protein kinase C epsilon type (PRKCE)

<b>Product Code</b>	CSB-EP018704HU-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q02156
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVVFNGLLKIKICEAVSLKPTAWSLRHAVGPRPQTFLLDPYIALNVDDSRIGQT ATKQKTNSPAWHDEFVTDVVCNGRKIELAVFHDAPIGYDDFVANCTIQFEELLQ NGSRHFEDWIDLEPEGRVYVIIDLSGSSGEAPKDNEERVFRERMRPRKRQGA VRRRVHQVNGHKFMATYLRQPTYCSHCRDFIWGVIGKQGYQCQVCTCVVHK RCHELIITKCAGLKKQETPDQVGSQRFSVNMPHKFGIHNYKVPTFCDHCGSLL WGLLRQGLQCKVCKMNVHRR CETNVAPNCGVDARGIAKVLADLGVTPDKITN SGQRRKLIAGAESPQPASGSSPSEEDRSKSAPTSPCDQEIKELENNIRKALS F DNRGEEHRAASSPDGQLMSPGENGEVRQGQAKRLGLDEFNFIKVLGKGSFG KVMLAELKGKDEVYAVKVLKQDVILQDDDDV DCTMTEKRILALARKHPYLTQLYC CFQTKDRLFFVMEYVNGGDLMFQIQRSRKFD EPRSRFYAAEVTSALMFLHQH GVIYRDLKLDNILLDAEGHCKLADFGMCKEGILNGVTTTTFCGTPDYIAPEILQE LEYGPSVDWWALGVLMYEMMAGQPPFEADNEDDLFESILHDDVLYPVWLSK EAVSILKAFMTKNPHKRLGCVASQNGEDA IKQHPFFKEIDWVLEQQKIKPPFK PRIKTKRDVNNFDQDFTREEPVLT LVDEAIVKQINQEEFKGFSYFGEDLMP
<b>Source</b>	E.coli
<b>Target Names</b>	PRKCE
<b>Expression Region</b>	1-737
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.