



Recombinant Human Lamin-B2 (LMNB2)

Product Code	CSB-BP013005HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q03252
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSPPPSPGRRREQRRPRAAATMATPLPGRAGGPATPLSPTRLRLQEKEELRE LNDRLAHYIDRVRALELENDRLLLKISEKEEVTTREVSIGIKALYESELADARRVL DETARERARLQIEIGKLRAELDEVNKSARKKREGELTVAQGRVKDLESLFHRSE VELAAALSDKRGLESDVAELRAQLAKAEDGHAVAKKQLEKETLMRVDLENRC QSLQEELDFRKSVEEEEVRETRRRHERRLVEVDSSRQQEYDFKMAQALEELR SQHDEQVRLYKLELEQTYQAKLDSAKLSSDQNDKAASAAREELKEARMRLES LSYQLSGLQKQASAAEDRIRELEEAMAGERDKFRKMLDAKEQEMTEMRDVM QQQLAEYQELLDVKLALDMEINAYRKLLEGEEERLKLSPSPSSRVTVSRATSS SSGSLSATGRLGRSKRKRLEVEEPLGSGPSVLGTGTGGSGGFHLAQQASAS GSVSIEEIDLEGKFVQLKNNSDKDQSLGNWRIKRQVLEGEEIAYKFTP KYILRA GQMVTVWAAGAGVAHSPSTLVWKGQSSWGTGESFRTVLVNADGEEVAMR TVKKSSVMRENGENGE EEEEEAEFG EEDLFHQGGDPRTT SRGC
Source	Baculovirus
Target Names	LMNB2
Expression Region	1-617
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.