



Recombinant Human Transducin-like enhancer protein 1 (TLE1)

Product Code	CSB-EP717264HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q04724
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFPQSRHPTPHQAAGQPFFKFTIPESLDRIKEEFQFLQAQYHSLKLECEKLASE KTEMQRHYVMYYEMSYGLNIEMHKQTEIAKRLNTICAQVIPFLSQEHQQQVAQ AVERAKQVTMAELNAIIGQQQLQAQHLSHGHGPPVPLTPHPISGLQPPGIPPLG GSAGLLALSSALSGQSHLAIKDDKKHHDAEHHRDREPGTNSNLLVPDSLRTGTD KRRNGPEFSNDIKRKVDDKDSSHYDSDGDKSDDNLVVDVSNEDPSSPRASP AHSPRENGIDKNRLLKKDASSSPASTASSASSTSLKSKEMSLHEKASTPVLKS STPTPRSDMPTPGTSATPGLRPLGKPPAIDPLVNQAAAGLRTPLAVPGPYPA PFGMVPHAGMNGELTSPGAAYASLHNMSQMSAAAAAAVVAYGRSPMVGF DPPPHMRVPTIPPNLGIPGGKPAYSFHVTADGQMMPVPFPDPALIGPIPRH ARQINTLNHGEVVCAVTISNPTRHVYTGKGKCVKVDISHPGNKSPVSQDCL NRDNYIRSKLLPDGCTLIVGGEASTLSIWDLAAPTPIKAELTSSAPACYALAI SPDSKVCFSCCSDGNIWDLHNQTLVRQFQGHDTGASCIDISNDGTLKWTG GLDNTVRSWDLREGRQLQQHDFTSQIFSLGYCPTGEWLAVGMESNNVEVLH VNKPDKYQLHLHESCVLSLKFAFCGKWFVSTGKDNLLNAWRTPYGASIFQSK ESSVLSCDISVDDKYIVTGSGDKKATVYEVY
Source	E.coli
Target Names	TLE1
Expression Region	1-770
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.