



Recombinant Human Glutamine--fructose-6-phosphate aminotransferase (GFPT1)

Product Code	CSB-BP009377HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q06210
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	CGIFAYLNYHVPRTTRREILETLIKGLQRLEYRGYDSAGVGFDDGGNDKDWEANA CKIQLIKKKGKVKALDEEVHKQQDMDLDIEFDVHLGIAHTRWATHGEPSPVNS HPQRSCKNNEFIVIHNGIITNYKDLKKFLESKGYDFESETDTETIAKLVKMYMDN RESQDTSFTTLVERVIQQLEGAFALVFKSVHFPGQAVGTRRGSPLLIGVRSEH KLSTDHIPILYRTARTQIGSKFTRWGSQGERGKDKKGCNLSRVDSTTCLFPV EEKAVEYYFASDASAVIEHTNRVIFLEDDDDVAAVVDGRLSIHRIKRTAGDHPGR AVQTLQMEQLQIMKGNFSSFMQKEIFEQPESVNTMRGRVNFDDYTVNLGGL KDHIKEIQRCRRLIACGTSYHAGVATRQVLEELTELPVMVELASDFLDRNTPV FRDDVCFFLSQSGETADTLMGLRYCKERGALTVGITNTVGSSISRETDCGVHI NAGPEIGVASTKAYTSQFVSLVMFALMMCDRISMQERRKEIMLGLKRLPDLIK EVLSDMDEIQKLATELYHQKSVLIMGRGYHYATCLEGALKIKEITYMHSEGILAG ELKHGPLALVDKLMPVIMIIMRDHTYAKCQNALQQVVARQGRPVICDKEDTET IKNTKRTIKVPHSVDCCLQGILSVIPLQLLAFHLAVLRGYDVFPRNLAKSVTVE
Source	Baculovirus
Target Names	GFPT1
Expression Region	2-699
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.